

DispensMate Bottle-Top Dispensers

Quick and easy volume adjustment

Smooth piston design provides effortless dispensing

Adapters allow fitting to standard reagent bottles

Closure cap protects, against reagent contact and is easy to attach and remove, even while wearing gloves

Filling tube adjusts easily to different size bottles

Features

- Excellent chemical resistance
- Fully autoclavable at 121°C
- Four ranges of bottle-top dispenser cover a volume range from 0.5mL to 50mL
- Easy to clean and maintain
- The optional flexible discharge tube with safety handle permits fast and precise dispensing
- Made of PTFE, FEP, BSG, PP
- Vapor pressure Max. 500mbar, viscosity max. 500mm²/s, temperature max. 40°C, density max. 2.2g/cm³
- DispensMate plus is supplied with S40, GL32, GL38, GL25, GL28 sized adapters

Specifications

Volume Range mL	Graduation mL	A _±		CV _≤	
		%	μL	%	μL
0.5-5	0.1	0.5	25	0.1	5
1.0-10	0.2	0.5	50	0.1	10
2.5-25	0.5	0.5	125	0.1	25
5.0-50	1.0	0.5	250	0.1	50

Accessories

S40, Adapter, 45/40mm

GL32, Adapter, 45/32mm

GL38, Adapter, 45/38mm

GL25, Adapter, 32/25mm

GL28, Adapter, 32/28mm

Reagent Bottle (Brown, 1L, Options)

Filling tube(m)

Discharge tube(m)