


Fixed Top LED Orbital Shaker



Features:

- * LED display for time and speed.
- * 320*265mm platform available for different vessels: petri dish, flask, etc.
- * Compact design, ideal for use inside incubator.
- * Provide gentle and efficient mixing

Technical Parameters:

Model	SK-O180-F
Orbital Diameter	20mm
Max. Load Capacity	3kg
Platform Size	320*265mm
Motor Type	DC motor
Motor Input Power	16W
Motor Output Power	10W
Speed Range	70-200rpm
Time Setting Range	1-1199min
Operation Mode	Continuous or timing
Ambient Temp. & Humidity	5-40°C, 80%RH
Protection Type	IP21
Power Supply	AC220V±10%, 50/60Hz
External Size (W*D*H)	330*270*110mm
Package Size (W*D*H)	410*380*260mm
Gross Weight	5.1kg