

Features

- Pipette controller for use with pipettes from 0.1 - 100mL
- Eight speeds are available for aspirate and dispense liquid
- Large LCD display provides visual confirmation of remaining battery charge and speed settings
- Easy for one-handed operation
- Lightweight, ergonomic designs allow longer, fatigue-free pipetting
- Efficient lithium-ion battery offers long runtime on each charge
- Low battery alert
- Powerful motor fills a 25mL pipette in <5 seconds
- Supplied with integral but replaceable 0.45µm filter
- Can be used while charging

Large LCD display provides the remaining battery life and speed settings



Select the desired coarse speed by press the button +/-

Comfortable aspirate and dispense triggers control fine volume speed

The dispense button is pressed only up to a point where a check can be felt, the dispensing function will be based on the force of gravity and the liquid will simply flow out

The 0.45µm integral filter and adapter are replaceable

Silicone adapter suitable for 0.1-100mL pipettes that can be fully autoclaved at 121°C



Levo Plus Pipette Filler



Specifications

Aspirate Speeds	8
Dispense Speeds	8 Gravity Dispense
Battery	Lithium-Ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2~3 Hours
Pipette Types	Glass or Plastic Pipette(0.1-100mL), Pasteur Pipettes
Filter	0.45 µm Hydrophobic



Two-position stand

